

Protein digestion by S-Trap™ micro MS sample preparation

Adapted from PROTIFI company

Introduction

The following protocol describes a method for generation of tryptic peptides from crude lysates for LC-MS analysis. The method allows analysis of detergent lysed cells and tissues. Therefore, it is particularly suitable for studying entire proteomes and fractions containing biological membranes.

1. Materials

Solutions and Reagents (use MS-grade water for all solutions)

- 1 M TEAB stock.
- 1X lysis buffer: 5% sodium dodecyl sulfate (SDS), 50 mM triethylammonium bicarbonate (TEAB), pH 8.5.
- Reductant: 120 mM Tris(2-carboxyethyl)phosphine (TCEP) in water.
- Alkylator: 500mM of chloroacetamide (CIAA) in water.
- Acidifier: Phosphoric acid diluted to 27.5% with water.
- Binding/wash buffer: 100 mM TEAB (final) in 90% methanol.
- Trypsin Stock: 1 µg/µl in 50 mM TEAB.
- Digestion buffer: 50 mM TEAB containing sufficient trypsin stock solution yield 1 µg trypsin per 10 µg sample (weight:weight) (or 5 ug maximum of trypsin for an IP digestion) in 20 µl 50 mM TEAB; no less than 1 µg trypsin per micro column.
- Elution buffer 1: 50 mM TEAB in water.
- Elution buffer 2: 0.2% formic acid in water.
- Elution buffer 3: 50% acetonitrile in water.

Equipment

1.5 ml or 2 ml lobind tubes.

2. Methods

2.1 Sample preparation

Lysis and thiol-reduction

There is no universal protocol to lyse any kind biological material but heating in the presence of high concentrations of SDS and reduction of disulfide bridges usually results in total or at least high extent of lysis.

Examples for total lysate preparation:

50 μ l HeLa cell pellet (one \emptyset 15 cm dish, about 2 mg of total proteins), are lysed by mixing with in 500 μ l of 1X lysis buffer and incubation at 95°C for 3 min. The DNA has to be sheared by sonication to reduce the viscosity of the sample. Before starting sample processing the lysate has to be clarified by centrifugation at 16,000 x g for 5 min. Quantify your protein extract using BCA kit.

2.2 Sample processing

1. In a lobind tube, mix up to 100 μ g of your protein extract with 1X lysis buffer to reach a volume of 92 μ l.
2. Add 4 μ l of reductant (final concentration 5 mM TCEP) and incubate at 55 °C for 15 min.
3. Add 4 μ l of alkylator (final concentration 20 mM CIAA) and incubate at RT for 10 min.
4. Add 15 μ l of acidifier, micros, to the 100 μ l sample (final concentration ~2.5% phosphoric acid). Vortex. This step is essential to completely denature proteins and trap them efficiently. The pH will be \leq 1. If the sample pH is not \leq 1, add additional phosphoric acid to reach pH \leq 1.
5. Add 660 μ l of binding/wash buffer to the sample and mix. Sufficient proteins may make the colloidal protein particulates visible, giving a translucent appearance. All sample – including any resulting colloid – must be transferred to the S-Trap. **Do not centrifuge at this point.**
6. Place an S-Trap micro-column in a 2 ml lobind tube for waste flow through. Apply sample to S-Trap column (180 μ l maximum per shot, and centrifuge at 4,000 g for 1 min per shot.). Transfer all sample including anything insoluble into the S-Trap. No column preequilibration is necessary.

7. Add 150 µl binding/wash buffer; centrifuge at 4,000 g for 1 min. Repeat 3 times and discard flow through as necessary. For best results, rotate the S-Trap micro units (like a screw or knob) 180 degrees between the centrifugations of binding and wash steps. This is especially important when using a fixed-angle rotor because the spin column does not experience homogenous flow. A mark on the outside edge during centrifugation makes it easy to track rotations. Additional wash(es) may be performed and should be performed if contamination is observed by mass spectrometry. You cannot over wash mammalian proteins with binding/wash buffer. If the sample has particularly hydrophobic contaminants, such as lipids in brain, bone marrow or adipose tissue, or wax in FFPE, bind protein, wash three times with 50% CHCl₃ /50% MeOH, filling the entire column each time, then perform three washes with standard methanolic binding/wash buffer as normal. After each centrifugation step, make sure that all added solution has gone through the S-Trap column.

8. Centrifuge S-Trap column at 4,000 g for 1 min to fully remove binding/wash buffer. Failure to do so may result in dripping.

9. Transfer S-Trap micro column to a clean 1.5 ml lobind tube for the digestion.

10. Add 20 µl of digestion buffer containing protease, most commonly trypsin or a trypsin/lys-C mix, at a 1:10 weight to weight (wt:wt) ratio into the top of the S-Trap. By example, for 50 µg of protein to digest, apply 5 µg of trypsin dissolved in 20 µl. Do not apply less than 1 µg of trypsin for effective digestion. 50 mM TEAB, TRIS or ammonium bicarbonate are recommended as digestion buffers for trypsin or trypsin/lys-C mixes. Visually confirm no air bubbles are present at the top the trap. Bubbles prevent the digestion buffer from entering the trap. If bubbles are present, flick the tube to remove them and/or centrifuge extremely briefly on a bench top centrifuge, returning any flow through to the top of the spin column. The S-Trap is hydrophilic and will absorb the digestion buffer. Careful observation shows that the applied digestion buffer “sinks” by a fraction of a mm when first applied. Do not damage the matrix with pipette tips. Other proteases than trypsin can be used, however they may require different digestion buffers, temperatures, pHes, metal cofactors, and/or times. S-Traps have been found to be compatible with Tryp-N, Lys-N, Lys-C, Arg-C, Glu-C, chymotrypsin, elastase and pepsin. Trypsin/lys-C mixes generally perform better than trypsin alone. S-Traps can also be used for glycomics with PNGase F. Mass spec compatible detergents such as Rapigest™ are compatible with S-Trap sample digestion and for some sample types have been observed to aid in digestion. ProteaseMAX™ should not be used at elevated digestion temperatures due to accelerated autolysis. Trypsin/lys-C mixes generally perform better than trypsin alone.

11. Cap the S-Trap **loosely** to limit evaporative loss. **Do not make an airtight seal with the cap.** An airtight seal will force the solution out of the S-Trap during incubation. Add parafilm to prevent evaporation.

12. Incubate overnight at 37 °C. Preferably use a water bath or stationary thermomixer. Ensure the entire column is exposed to heat. Do not shake.

13. Add 40 µl of elution buffer 1 to the S-Trap then centrifuge (4,000 g, 1 min). Do not centrifuge the digestion through before applying elution buffer 1. Apply elution buffer 1 directly into the trap containing the digestion buffer that was incubated.

14. Add 40 µl of elution buffer 2 to the S-Trap then centrifuge (4,000 g, 1 min).

15. Add 40 μ l of elution buffer 3 to the S-Trap then centrifuge (4,000 g, 1 min). This elution assists in recovery of hydrophobic peptides. Other organics may also be used as needed.

16. Pool eluted peptides and speedvac to dryness.

17. Resuspend in 100 μ l 1% formic acid in MS-grade water. No further C18 cleanup is necessary.

2.3 Yield determination

Concentration of the peptides can be estimated by Nanodrop absorbance at 205 nm.