

# ImmunoPrecipitation using GFP-TRAP

Boisvert Lab, December 2020

- Resuspend cell pellet in cold lysis buffer (1 ml lysis buffer/100 mm petri dish).
- Allow lysis to continue 10 minutes on ice.
- Spin down insoluble material 10 minutes at 13,000 rpm, or until no pellet is detected after centrifugation.
- Add GFP-TRAP sepharose to lysate (10-20  $\mu$ l of washed beads per 100 mm petri dish).
- Tumble 2 hrs at 4°C.
- Spin down resin at 2,000 rpm for 3 minutes.
- Wash three times with lysis buffer.
- Wash twice with PBS 1X.
- Spin down resin and aspirate all liquid.
- Add Laemmli loading buffer to elute.
- Boil sample for 5 minutes.
- Load on SDS-PAGE.

## **Solutions:**

### Lysis Buffer

For 50 ml

1% Triton X-100  
150 mM NaCl  
20 mM Tris-HCl, pH 7.5  
0.1 mg/ml PMSF

5 ml of 10%  
1.5 ml of 5 M  
1 ml of 1 M  
500  $\mu$ l of 10  $\mu$ g/ $\mu$ l